



Benefit™

Triple Head Toothbrush

Brushes teeth with little or no wrist movement and gets better results!

Brushes Right Every Time!



Patented Triple Head Toothbrush is easier to keep on the teeth and enables the brushing of all sides of each tooth at once.

Available In Two Styles

- Benefit Toothbrush - standard handle for use with or without a cuff.
- Benefit Plus Toothbrush - with a wide, padded, easy to hold handle.

Both styles are great for people with arthritis, Parkinson's, MS and other problems. The padded handle provides the best grip.

- Benefit toothbrushes get the job done faster, easier and better.
- 3 flexible heads widen and contract to maintain the proper pressure on all surfaces of the tooth with every stroke.
- Side brushes are at an optimal 45 degree angle, which reduces the need to twist the arm and wrist. Softer outer bristles massage the gums, helping to prevent periodontal disease.
- Made in the USA.
- US Patents 5,305,491; 5,407,254; 5,673,545, D372,586.



Benedent Corporation 7795 Locust Place, Port Tobacco, Maryland 20677 ©Benedent Corp.

Tel: (800) 450-4977 & (301) 934-1204 Fax: (301) 934-0847 e-mail: rogerbayer@benedent.com www.benedent.com